

Product Datasheet

Anti-Phospho-TrkA(Tyr496) +TrkB(Tyr516)/Cy5.5 (orb1047206)

Description	Anti-Phospho-TrkA(Tyr496) +TrkB(Tyr516)/Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Cy5.5
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TrkA around the phosphorylation site of Tyr496
Target	Phospho-TrkA (Tyr496) +TrkB(Tyr516)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04629
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.