

Product Datasheet

CTGF antibody (orb10474)

Description

CTGF antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CTGF. Please contact us for the exact immunogen sequence. The peptide is available as orb374719.

Target

CTGF

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW


36 kDa

Uniprot ID
[Q5M8T4](#), [P29268](#)
NCBI
[21](#), [001901](#), [001892](#), [11](#)
Entrez
[1490](#)
Dilution Range


WB: 1:100-500, IHC-P: 1:100-800 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P staining of rat brain tissue using...



Immunohistochemical staining of paraffin...



IHC-P staining of rat trachea tissue usi...