



Product Datasheet

CTLA4 antibody (orb10475)


Description	CTLA4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CTLA-4 (141-223/223 aa)
Target	CTLA4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	P16410
Dilution Range	WB=1:200-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver



Immunohistochemical staining of Mouse Colon



Immunohistochemical staining of Mouse Tumor