

Product Datasheet

**Anti-phospho-PI3 kinase p85 alpha + gamma
(Tyr467 + Tyr199)/Cy5.5 (orb1047624)**

Description

Anti-phospho-PI3 kinase p85 alpha

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Cy5.5

Tested

FC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467

Target

phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
P27986
Dilution Range

Flow-Cyt=1ug/test,IF=1:100-500

Expiration Date

12 months from date of receipt.