



Product Datasheet

Connexin 32 antibody (orb10481)


Description	Connexin 32 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Connexin-32 (201-283/283 aa)
Target	Connexin-32
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	P08034
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot of rat cardiac muscle tissue...



IHC-P of rat cardiac muscle tissue (orb1...