

Product Datasheet

SARS-CoV-2 (2019-nCoV) Spike Neutralizing antibody (orb1048254)

Description	SARS-CoV-2 (2019-nCoV) Spike Neutralizing antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Spike S1 Protein (14-685/1213 aa)
Target	SARS-CoV-2 (2019-nCoV) Spike Neutralizing
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	Lot Dependent
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	140 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.