

Product Datasheet

Tau isoform 4 antibody (orb1048341)

Description	Tau isoform 4 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Tau isoform 4 protein
Target	Tau isoform 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	37/79 kDa
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.