



Product Datasheet

GP1BA antibody (orb1048342)

Description	GP1BA antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GP1BA/CD42b
Target	GP1BA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	67 kDa
Uniprot ID	P07359
Dilution Range	WB=1:500-2000, IHC-P=1:200-1000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse
spleen using
CD42b...



IHC-P
analysis of
rat spleen
using
CD42b...