

Product Datasheet

Anti-Phospho-FAK (Tyr407)/RBITC (orb1048627)

Description	Anti-Phospho-FAK (Tyr407)/RBITC
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	RBITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FAK around the phosphorylation site of Tyr407
Target	Phospho-FAK (Tyr407)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Uniprot ID	Q05397
Dilution Range	Flow-Cyt=1µg /Test,IF=1:100-500
Expiration Date	12 months from date of receipt.