



Product Datasheet

CXCR1 antibody (orb10487)


Description	CXCR1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CXCR1. Please contact us for the exact immunogen sequence. The peptide is available as orb12890.
Target	CXCR1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	40 kDa
Uniprot ID	Q810W6 , P25024 , P70612
NCBI	000625
Entrez	3577
Dilution Range	WB: 1:100-800, IHC-P: 1:100-600, IF/ICC: 1:100-600
Expiration Date	12 months from date of receipt.



IF image of guinea pig skin tissue using...



IHC-P image of rat lymph node tissue using...



Western blot analysis of rat spleen (lanes 1-6) using CXCR1 antibody.