


Product Datasheet

PCDHGB6 antibody (orb1049)

Description	PCDHGB6 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PCDHGB6 (21-120/930 aa)
Target	PCDHGB6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	98 kDa
Uniprot ID	Q9Y5F9
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Rsc96 cell
u...