

Product Datasheet

Anti-Phospho-Histone H3 (Ser10)/Cy5.5 (orb1049589)

Description	Anti-Phospho-Histone H3 (Ser10)/Cy5.5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Cy5.5
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Ser10
Target	Phospho-Histone H3 (Ser10)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	Flow-Cyt=0.2µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.