

Product Datasheet

Anti-phospho-Erk1 (Thr202 + Tyr204) /Cy7 (orb1050919)

Description	Anti-phospho-Erk1 (Thr202 + Tyr204) /Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy7
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat ERK1 around the phosphorylation site of Thr201/204
Target	phospho-Erk1 (Thr202 + Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.