



Product Datasheet

CK-HMW antibody (orb10510)

Description	CK-HMW antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Cytokeratin-13 (251-350/458 aa)
Target	Cytokeratin 13
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	P13646
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Tongue(Mouse);Primary:Anti-
...



Western blot:HeLa (Human)
Cell;Primary:A...