

Product Datasheet

Anti-phospho-JNK1 + 2 + 3 (Thr183+Tyr185)/RBITC (orb1051135)

Description	Anti-phospho-JNK1 + 2 + 3 (Thr183+Tyr185)/RBITC
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	RBITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JNK1 around the phosphorylation site of Thr183/Tyr185
Target	phospho-JNK1 + 2 + 3 (Thr183+Tyr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45983
Dilution Range	Flow-Cyt=1µg /test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.