

Product Datasheet

Anti-phospho-ERK1 + 2 (Thr183/Tyr185)/Cy7 (orb1051462)

Description	Anti-phospho-ERK1 + 2 (Thr183/Tyr185)/Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Cy7
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse ERK1 around the phosphorylation site of Thr183/Tyr185
Target	phospho-ERK1 + 2 (Thr183/Tyr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.