

Product Datasheet

Dopamine D2 receptor antibody (orb10515)

Description

Dopamine D2 receptor antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Dopamine D2 receptor. Please contact us for the exact immunogen sequence. The peptide is available as orb12891.

Target

Dopamine D2 receptor

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

51 kDa

Uniprot ID

[P61168](#), [P14416](#), [P61169](#)

NCBI

[31](#), [000786](#), [000795](#), [11](#)

Entrez


[1813](#)

Dilution Range


WB: 1:100-2000, IHC-P: 1:100-500, IF/ICC: 1:100-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunofluorescence analysis of rat brain...



IHC-P staining of rat brain tissue using...