



Product Datasheet

Connexin 26 antibody (orb10516)


Description	Connexin 26 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Connexin-26 (81-180/226 aa)
Target	Connexin 26
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26 kDa
Uniprot ID	Q00977
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot of rat ovary tissue lysate ...



IHC-P of rat mouse ovary tissue (orb1015...