

Product Datasheet

Anti-BrdU/Cy5.5 (orb1052421)

Description	Anti-BrdU/Cy5.5
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Cy5.5
Tested Applications	IF
Immunogen	KLH conjugated Brdu
Target	BrdU(Proliferation Marker)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
CAS Number	59-14-3
Expiration Date	12 months from date of receipt.