

Product Datasheet

Anti-phospho-GSK-3 Beta(Ser9)/Cy7 (orb1052879)

Description	Anti-phospho-GSK-3 Beta(Ser9)/Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy7
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Ser9
Target	phospho-GSK-3 Beta (Ser9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49841
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.