

Product Datasheet

Anti-Substance P/RBITC (orb1053276)

Description	Anti-Substance P/RBITC
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	RBITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Substance P (RPKPQQFFGLM)
Target	Substance P
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20366
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.