

Product Datasheet

Anti-HMGB1/Cy7 (orb1053955)

Description	Anti-HMGB1/Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy7
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human HMGB1
Target	HMGB1
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09429
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.