



Product Datasheet

Dnmt1 antibody (orb10546)


Description	Dnmt1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Dnmt1 (21-120/1616 aa)
Target	Dnmt1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	178 kDa
Uniprot ID	P26358
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical analysis of formalin...



IHC-P of human lung carcinoma (Dnmt1 ant...