



Product Datasheet

DNMT3B antibody (orb10548)

Description	DNMT3B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Dnmt3 Beta (1-80/853 aa)
Target	Dnmt3b
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	94 kDa
Uniprot ID	Q9UBC3
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HepG2 cell u...



Immunohistochemical analysis of Mouse st...