

Product Datasheet

Anti-Estrogen Receptor alpha + beta/Cy7 (orb1054854)

Description	Anti-Estrogen Receptor alpha + beta/Cy7
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Cy7
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Estradiol Receptor alpha + beta
Target	Estrogen Receptor alpha + beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.