


Product Datasheet

DR4 antibody (orb10553)

Description	DR4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human TNFRSF10A protein
Target	TNFRSF10A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	000220
Dilution Range	WB=1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...