

Product Datasheet

**Anti-Phospho-RPS6 (Ser240 + Ser244)/RBITC
(orb1055489)**

Description

Rabbit polyclonal antibody to Anti-

Ser244)/RBITC

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

RBITC

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human S6 Ribosomal Protein around the phosphorylation site of Ser240/244

Target

RPS6

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

29kDa

Uniprot ID
P62753
Expiration Date

12 months from date of receipt.