


Product Datasheet

MMP9 antibody (orb10565)

Description	Rabbit polyclonal antibody to MMP9
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human MMP9(NP_004985.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	92KDa
Uniprot ID	P14780
Dilution Range	WB: 1:500 - 1:2000, IHC-P: 1:50 - 1:200, IF/ICC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...