


Product Datasheet

HPV16 E7 antibody (orb10573)


Description	HPV16 E7 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from HPV16 E7 protein (41-90/98 aa)
Target	HPV16 E7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P03129
Dilution Range	Flow-Cyt=1ug/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human uterus tissue (HPV16 E7 a...



Immunohistochemical analysis of formalin...