

Product Datasheet

Anti-Insulin receptor subunit beta/Cy5.5 (orb1057425)

Description	Anti-Insulin receptor subunit beta/Cy5.5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Insulin receptor subunit beta
Target	Insulin receptor subunit beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.