



Product Datasheet

EGFR antibody (orb10580)


Description	EGFR antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human EGFR (951-1050/1210 aa)
Target	EGFR
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130 kDa
Uniprot ID	P00533
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P of human gastric cancer tissue (EG...)



Immunohistochemical analysis of paraffin...



Paraffin-embedded human lung carcinoma t...