

Product Datasheet

Anti-phospho-beta Actin (Tyr53)/RBITC (orb1058016)

Description	Anti-phospho-beta Actin (Tyr53)/RBITC
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	RBITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human beta Actin around the phosphorylation site of Tyr53
Target	phospho-beta Actin (Tyr53)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60709
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.