



Product Datasheet

Emp1 antibody (orb10588)


Description	Emp1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse Emp1.Å
Target	Emp1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	18 kDa
Uniprot ID	P54848 , P47801
NCBI	034258 , 010128
Entrez	13730
Dilution Range	WB: 1:100-800, IHC-P: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



IHC-P staining of rat liver tissue using...



WB analysis of rat lung (lane 1), rat li...



IHC-P image of human endometrial cancer ...