

Product Datasheet

HSP20 (Phospho-Ser16) antibody (orb105917)

Description	Rabbit polyclonal to HSP20 (Phospho-Ser16). HSP20 a small heat shock protein of the HSP20 family. No...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...