


Product Datasheet

Chk1 (Phospho-Ser286) antibody (orb105919)

Description	Rabbit polyclonal to Chk1 (Phospho-Ser286). DNA damage induced protein phosphorylation; regulation O...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...