


Product Datasheet

MAP3k7 (Phospho-Thr187) antibody (orb105931)

Description	Rabbit polyclonal to MAP3k7 (Phospho-Thr187). AK1 a protein kinase of the MLK family. Mediates signa...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...