


Product Datasheet

JAk2 (Phospho-Tyr221) antibody (orb105935)

Description	Rabbit polyclonal to JAk2 (Phospho-Tyr221). This gene product is a protein tyrosine kinase involved ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...