


Product Datasheet

Dab1 (Phospho-Tyr232) antibody (orb105939)

Description	Rabbit polyclonal to Dab1 (Phospho-Tyr232). DAB1 an adaptor molecule functioning in neural developme...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of LOVO whole cell...