


Product Datasheet

MSK1 (Phospho-Ser376) antibody (orb105940)

Description	Rabbit polyclonal to MSK1 (Phospho-Ser376). an AGC kinase of the RSK family. Has two kinase domains ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...