


Product Datasheet

SNAI1 (Phospho-Ser246) antibody (orb105942)

Description	Rabbit polyclonal to SNAI1 (Phospho-Ser246). This protein has many roles during postimplantation dev...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HT29 whole cell...