


Product Datasheet

Chk2 (Phospho-Ser516) antibody (orb105945)

Description	Rabbit polyclonal to Chk2 (Phospho-Ser516). In response to DNA damage and replication blocks, cell C...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...