


Product Datasheet

EGFR (Phospho-Thr693) antibody (orb105951)

Description	Rabbit polyclonal to EGFR (Phospho-Thr693). EGFR is a receptor tyrosine kinase. Receptor for epiderm...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 whole cell...