


Product Datasheet

Estrogen Receptor-alpha (Phospho-Ser106) antibody (orb105965)

Description	Rabbit polyclonal to Estrogen Receptor-alpha (Phospho-Ser106). ER-alpha is a nuclear hormone recepto...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 whole cell...