


## Product Datasheet

# Glucocorticoid Receptor (Phospho-Ser203) antibody (orb105971)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Glucocorticoid Receptor (Phospho-Ser203) |
| <b>Species/Host</b>        | Rabbit                                                                 |
| <b>Reactivity</b>          | Human, Mouse                                                           |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | WB                                                                     |
| <b>Note</b>                | For research use only                                                  |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Dilution Range</b>      | WB:1:500-1:2000                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Western blot analysis of A549 whole cell...