


Product Datasheet

p53 (Phospho-Ser15) antibody (orb105976)

Description	Rabbit polyclonal antibody to p53 (Phospho-Ser15)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...