


Product Datasheet

HSP27 (Phospho-Ser82) antibody (orb105982)

Description	Rabbit polyclonal to HSP27 (Phospho-Ser82). HSP27 is a small heat shock protein that is regulated bo...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...