


Product Datasheet

Integrin beta3 (Phospho-Tyr785) antibody (orb105985)

Description	Rabbit polyclonal to Integrin beta3 (Phospho-Tyr785). ITGB3 integrin alpha-V/beta-3 is a receptor fo...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...