


## Product Datasheet

# Progesterone Receptor (Phospho-Ser400) antibody (orb106003)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to Progesterone Receptor (Phospho-Ser400). Progesterone plays a central role in th... |
| <b>Species/Host</b>        | Rabbit                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | ICC, IF, WB                                                                                             |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                              |
| <b>Dilution Range</b>      | WB:1:500-1:2000, IF/ICC:1:100-1:500                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Western blot analysis of 293 whole cell ...