

Product Datasheet

EPOR antibody (orb10601)

Description

EPOR antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC-P, WB

Immunogen

KLH conjugated synthetic peptide derived from human EPOR. Please contact us for the exact immunogen sequence. The peptide is available as orb12893.

Target

EPOR

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

55 kDa

Uniprot ID

[P19235](#), [P14753](#), [Q07303](#)

NCBI

[000112](#), [31](#), [000121](#), [11](#)

Entrez


[2057](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:200-1000, IF/ICC: 1:100-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunofluorescence analysis of human lun...



Immunohistochemical staining of human lu...