

Product Datasheet

Histone H2A (Phospho-Thr121) antibody (orb106024)

Description	Rabbit polyclonal to Histone H2A (Phospho-Thr121). H2A.1 Core component of nucleosome. Nucleosomes w...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.